


smooth[®]

TECHNICAL INFO INFORMAZIONI TECNICHE



MAX LOAD 120 Kg

Attention it could change the setting of the wheelchair.



Before mounting, read carefully the technical info and verify the dimensions.

For a correct assembly, it is recommended to mount the fork in order to obtain the angle ① and the angle ② at 90°

CARICO MASSIMO 120 Kg

Attenzione potrebbe cambiare l'assetto della carrozzina




Prima del montaggio, leggere attentamente le informazioni tecniche e verificare le dimensioni.

Per un corretto montaggio si consiglia di montare la forcella al fine di ottenere l'angolo ① e l'angolo ② a 90°

POSSIBLE CONFIGURATIONS POSSIBILI CONFIGURAZIONI

The forks Smooth, can be mounted with 4 or 5 inch wheels

Le forcelle Smooth, possono essere montate con ruote da 4 o 5 pollici



ARIA WHEELS SRLS

Registered Office_Via della Divisione Torino, 92 00143 Roma - Italy

Production Dept._Via A.Volta, 7/G 36030 Costabissara (VI) - Italy


VAT IT 12848301003 | tel.+39 04441322533

www.ariawheels.com | info@ariawheels.com


s//ooth[®]

COMPATIBILITY WITH THE ARIA[®] PRODUCTS COMPATIBILITA' CON I PRODOTTI ARIA[®]

STANDARD FRAME TELAIO STANDARD




SHF480 max KHL400
SHF490 max KHL410
SHF500 max KHL420
SHF510 max KHL430
SHF520 max KHL440




SHF490 max KHL380
SHF500 max KHL390
SHF510 max KHL400
SHF520 max KHL410
SHF530 max KHL420

LOWERED FRAME TELAIO BASSATO



SHF450 max KHL360
SHF460 max KHL370
SHF470 max KHL380
SHF480 max KHL390
SHF490 max KHL400



SHF460 max KHL340
SHF470 max KHL350
SHF480 max KHL360
SHF490 max KHL370
SHF500 max KHL380

ARIA WHEELS SRLS

Registered Office_Via della Divisione Torino, 92 00143 Roma - Italy
Production Dept._Via A.Volta, 7/G 36030 Costabissara (VI) - Italy
VAT IT 12848301003 | tel. +39 04441322533
www.ariawheels.com | info@ariawheels.com